

MICRO-SUCTION

Wax Removal Equipment List

For pricing email – sales@jrmedical.co.nz

O SCOPE



SPIGGLE & THEIS

Medizintechnik

VOROTEK

experience. knowledge. performance








essential ear microsuction training









JR Medical







Code	Description	Qty	Image
VSOH-3D	<p>VorOtek O Scope 3 Diopter</p> <p>Head Band or Spec Frame</p> <p>12mm Converged Binocular Optics & Fully Integrated LED Illumination with 3 dioptre Magnification. Adjustable IPD. Includes Power Pack, Charger, Carry Case,</p>	1	
VSOH-3D-Kit	<p>VorOtek O Scope - HeadBand or Spec Frame - 3D WITH STARTER KIT</p> <p>12mm Converged Binocular Optics & Fully Integrated LED Illumination with 3 dioptre Magnification. Adjustable IPD. Includes Power Pack, Charger, Carry Case,</p> <ul style="list-style-type: none"> - 2 x Curettes 2.5 & 3.5mm - 2 x 90 hooks S & L - Wool Carrier - 15 x speculums XS,S,M,L,XL - Pneumatic Oto scope 6 x diopter with puffer 	1	
VEBL-10	<p>Vorotek Battery Pack 10 Hours</p>	1	






Code	Description	Qty	Image
VECM-10	Vorotek Charger for Battery Pack	1	
LTA270	Aspimed Surgical Pump 2.3 Piston Pump Power supply 230V/50Hz Power consumption 110VA Continuous use 1000cc jar Maximum suction flow 40l/min Maximum suction -0,80 bar High vacuum – High flow Includes all required accessories to get started	1	
LRA372	Replacement jar complete with lid with 1000 cc anti-overflow float, transparent for Aspimed 2.3	1	
LRA204	Replacement filter, antibacterial, transparent Aspimed 2.3	1	
LRA120	Replacement tube, 6 x 10 mm transparent silicone for Aspimed 2.3	1	




Code	Description	Qty	Image
DEF-1186	Suction Tubing Extender 5mm x 9mm x 30cm Shorter length tubing to assist easing the handling and weight of suction tube line Includes 3 x Wrapped Sterile Can re-use	1	
Suctions Tips 303025	Single use suction tips Spiggle & Theis STERILE Wrapped Box of 50 303025 – 2.5 x 80mm Blue	1	
Suction Tips 302020	Single use suction tips Spiggle & Theis STERILE Wrapped Box of 50 302020 – 2.0 x 70mm Green	1	
Suctions Tips 302016	Single use suction tips Spiggle & Theis STERILE Wrapped Box of 50 302016 – 1.6 x 70mm White	1	
Suctions Tips 302014	Single use suction tips Spiggle & Theis STERILE Wrapped Box of 50 302014 – 1.4 x 70mm Red	1	



Code	Description	Qty	Image
Suctions Tips 302012	Single use suction tips Spiggle & Theis STERILE Wrapped Box of 50 302012 – 1.2 x 70mm Pink	1	
Suctions Tips 302010	Single use suction tips Spiggle & Theis STERILE Wrapped Box of 50 302010 – 1.0 x 70mm Cream	1	
Suctions Tips 302007	Single use suction tips Spiggle & Theis STERILE Wrapped Box of 50 302007 – 0.7 x 70mm Black	1	
NZ323	Adept Suction Controllers Box of 20	1	
10-308-06-05 10-308-08-05 10-308-10-05 10-308-12-05 10-308-15-05 10-308-18-05 10-308-20-05 10-308-22-05 10-308-25-05	German Made Suction Tips Re-Usable - Autoclavable BUNDLE OF 5x Size Range 0.6mm 0.8mm 1.0mm 1.2mm 1.5mm 1.8mm 2.0mm 2.2mm 2.5mm		

Code	Description	Qty	Image
10-306-05	House suction handle with suction control 5 cm Autoclavable Re-usable German Made	1	
Wax Rings / Curettes	Vorotek Stainless Steel, Hand Made Sizes 2.0mm 2.5mm 3.0mm 3.5mm 4.0mm	1	
Wool Carrier	Vorotek Wool Carrier Stainless Steel	1	
90 Degree Hooks	Vorotek Surgical grade Stainless Steel. Operating room quality. Reusable Available in 2 sizes short or long	1	
Pneu Otoscope	Pneumatic Otoscope	1	
Vorotek Speculum	Vorotek Autoclavable Plastic Sizes XS, S, M, L, XL	1	

Code	Description	Qty	Image
Speculum Set	Vorotek Set of 15 x Vorotek Speculum 3 x XSmall 3 x Small 3 x Medium 3 x Large 3 x XLarge	1	
Speculum Wall Mount	VorOtek Speculum Dispenser (Wall mounted) Including 100 Speculums	1	
3.5 Crocs	Adept Medical Grade Low Stainless Steel Jaw 3.5mm, 75mm	1	
5001007	Single Use Disposable cerumen hooks 0.3mm, 13.5cm length clean – non sterile (box of 100x) Spiggle & Theis GmbH	1	
5001008	Single Use Disposable Double End Ear Curette with Cotton clean – non sterile (box of 100x) Spiggle & Theis GmbH	1	

Code	Description	Qty	Image
CliniMax Starter kit	<p>Instrument Detergent Includes</p> <ul style="list-style-type: none"> - 5 litres Clinimax - Clinimax 30ml Pump for 5L bottle - Clinimax 500ml Empty Bottle - Clinimax 4ml Pump for bottle - Clinimax White Spray Trigger <p>Mildly alkaline, non-hazardous, multipurpose detergent for Medical and Dental applications</p> <ul style="list-style-type: none"> • Non-hazardous • Odourless • Low foaming (to allow the user to see what is happening in solution) • Bacteriostatic (prevents the growth of bacteria) 		
10-276-02	<p>Hartmann ear forceps serrated, straight, fine 6.0x 0.8mm, 8 cm working length</p> <p>Spiggle & Theis GERMAN HIGH GRADE STAINLESS STEEL</p>		
10-276-00	<p>Hartmann ear forceps straight, serrated 8 x1.8 mm, length 8 cm</p> <p>Spiggle & Theis GERMAN HIGH GRADE STAINLESS STEEL</p>		



What will you learn?



How to develop Theoretical and Practical Skills using Best Practice Technique for Ear Microsuction and Curetage



How to Manage Problematic Cerumen and to Diagnose/Treat those Conditions and Diseases



How to understand when to refer certain Ear Conditions to other Professionals



The course is designed for General Practitioners; Registered or Enrolled Nurses and Audiologists

Meet your *trainer*

Brenda Lockie-Julian RN, BA, MBA

Registered Nurse, Vorotek Product Specialist, Ear Nurse Clinician Trainer

Brenda is an Australian trained ear nurse associated with Vorotek since 2010. She has been instrumental in setting up ear clinics in Auckland and has trained and mentored GPs, Registered Nurses and Audiologists throughout New Zealand for several years. Brenda is passionate about sharing her product and clinical knowledge and providing ongoing support for trainees.



Got Questions? Please Email me

brenda@eemt.co.nz