



## PENTAX Medical

### ENT Scopes

PENTAX Medical  
ENT endoscopes,  
with high-quality optics  
and a diverse range  
of features, help  
you identify tissue  
anomalies and  
improve detection  
of subtle pathology

**PENTAX**  
**MEDICAL**  
*Excellence in Focus*



## Rigid Laryngoscopes

The gold standard for stroboscopy, providing bright, crystal-clear images

- Integrated fiber bundle light guide for superior light transmission
- Patient comfort provided by gently tapered end
- Ergonomic design for transoral endoscopy



Model	Total length (mm)	Outer diameter (mm)	Working length (mm)	Actual Visual Field Angle (°)	Angled View (°)	Highlights
9106	252	10	189	35	70	Standard diameter
9108	225	6	189	35	70	Smaller diameter

## Fiber Nasolaryngoscopes

For cost-effective transnasal exams

- High-quality optics with excellent light transmission
- Extremely flexible for patient comfort
- Maximum mobility with the available portable LED light source for the FNL -10RBS



Model	Distal Tip (mm)	Channel (mm)	Working length (mm)	Angle of View (°)	Angulation Up/Down (°)	Highlights
FNL - 7RP3	2.4	N/A	300	75	130/130	Small diameter for pediatric use
FNL - 10RP3	3.4	N/A	300	75	130/130	Our best-selling fiber scope
FNL - 15RP3	4.8	2.2	300	75	130/130	Working channel
FNL - 10RBS	3.4	N/A	300	75	130/130	Portable LED light source

## Video Nasolaryngoscopes

The highest resolution scopes, including a high-definition option

- High-resolution CCD provides sharp, clear images
- Bifurcated light cable for use with stroboscopic light source
- Customizable buttons put commonly used settings at your fingertips



Model	Distal Tip (mm)	Channel (mm)	Working length (mm)	Angle of View (°)	Angulation Up/Down (°)	Highlights
VNL -1070STK	3.1	N/A	300	85	130/130	Small diameter for pediatric use
VNL -1190STK	4.1	N/A	300	80	130/130	Our best-selling video scope
VNL -1570STK	4.9	2.0	300	85	130/130	Working channel
VNL -1590STi	5.1	N/A	300	100	130/130	HD resolution
EE-1580K	5.5	2.0	600	140	210/120	Esophagoscope

*Note: All videoscopes require a PENTAX video processor. Please contact your product specialist for information on compatible processors.*

### PENTAX Medical

3 Paragon Drive

Montvale, NJ 07645

Phone +1 800 289 5297 • Fax +1 201 391 2063 • [www.kaypentax.com](http://www.kaypentax.com)

**PENTAX**  
**MEDICAL**

Copyright 2013. All Rights Reserved. All company and product names and marks contained within are federally registered trademarks, trademarks or service marks of PENTAX of America, Inc. Product specifications are subject to change without notice, and without any obligation on the part of the manufacturer.

MK-436 Rev: A